BAUMER

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Termodesinfectora
Modelo: TW-E2000

Código: MM.01.09.0006

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Folha: 1 de 75

Data: 28/11/2019

BAUMER

TERMODESINFECTORA

Registro ANVISA nº: 10345500051

Modelo: TW-E2000

M A N U A L D E M A N U T E N Ç Ã O

Código

Modelo

Fabricado em

Núm. Série

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Maintenance Manual

Washer-disinfector
Model: TW-E2000

Code: MM.01.09.0006
Version / Revision: 2.001
Page: 39 of 75
Date: 28/11/2019

BAUMER

WASHER-DESINFECTOR

ANVISA Registry Nº: 10345500051

Model: TW-E2000

M A I N T E N A N C E M A N U A L

Code:

Model:

Made in

Serial Number:

Baumer - Commitment to health
Technical responsible: Eng. Sérgio Yukio Koseki
CREA-SP: 0601577094 - Cart: 157709/D
Start of Technical Responsibility: 1/31/1994
Manufactured and Distributed by Baumer S.A.
STIC Division - Sterilization and Infection Control
Prefeito Antonio Tavares Leite Ave, 181 • Parque da Empresa
Mail Box 1081 • 13803-330 • Mogi Mirim • SP • Brazilian Industry
E-mail: baumer@baumer.com.br

Congratulations!

You have purchased a Baumer product and we are committed to matching your trust. To ensure the best performance of your product, read carefully and follow the instructions in this document.

We emphasize that the photos, figures and drawings are illustrative, being subject to variations without previous notification.

Baumer S.A. considers itself responsible for the safety, reliability and performance of its product provided that:

- The necessary supplies are available and meet the specifications of this guide;
- Only Baumer S.A. or authorized agents act in repairs and maintenance of the equipment;
- The product is used in accordance with the Instructions for Use, Installation Guide and Maintenance Manual supplied with the equipment.

Baumer S.A. is not liable for damages caused during the transport of your product. It is the responsibility of the buyer to make a conference at the time of delivery, registering any occurrence in the document of receipt of the product, immediately triggering the carrier and / or the insurer in case of damages, and informing Baumer S.A. of the event.

Baumer S.A., through its service channels, can monitor the receipt and withdrawal of its products from the original packaging. Consult us on how to hire these services.

The National Sanitary Regulation defines the need to perform the Installation Qualification (IQ) and Operational (OQ) of this equipment. These qualifications must be performed by the equipment supplier or by your service channels, after installation and repeated every year or after important interventions in the equipment.

Usually these services are negotiated with the buyer at the time of purchase of the equipment. If this has not happened, consult us on how to hire these services.

Also a Performance Qualification (PQ), to validate the processes performed by the equipment, with the parameters defined by the user, is a requirement of sanitary regulation. Baumer S.A. has availability to perform this service for its clients or to monitor the performance of this qualification when performed by third parties. Consult us on how to hire these services.

Baumer S.A. and its service channels have other services that can be contracted by users of our products, which guarantee an extended service life and a perfect functioning of this equipment, minimizing occurrences that may interrupt the use of the same. Consult us for contracting preventive maintenance services and specialized technical support.

We remain at your disposal for further clarification and we are sure that this product will be useful and productive for this institution for many years.

Baumer S.A. reserves the right to make, without prior notification, modifications to the product mentioned in this document, in compliance with the policy of continuous development and improvement of our products.

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without the prior permission of Baumer S.A.. For further clarification or contracting of these services, please contact MEDIXX, a company of Baumer S.A. group responsible for after-sales management of our equipment.

Summary

1. Introduction	42
2. Baumer Limited Warranty Term	43
3. Replacement of Safety Items.....	45
3.1. Replacing the HMI battery	45
3.2. Replacing the fuses.....	45
3.3. Replacing the emergency button	46
3.4. Replacing the general circuit breaker	46
3.5. Replacing the unipolar circuit breaker	47
3.6. Replacing the motor circuit breaker.....	48
3.7. Replacing the air pressure switch	48
3.8. Replacing the air pressure switch on the channel	49
3.9. Replacing the air pressure switch	49
3.10. Replacing the anti-crushing system.....	50
3.11. Replacing the thermostat the water resistances.....	51
3.12. Replacing the air resistances thermostat.....	51
4. Other maintenance.....	52
4.1. Door gasket replacement and maintenance	52
5. Preventive Maintenance Schedule	56
5.1. Periodic preventive maintenance on schedule:	56
5.1.1. Quarterly:	56
5.1.2. Semiannually:.....	56
5.1.3. Annually:	56
6. Spare parts	57
7. Calibration and Maintenance Routines	59
7.1. Manual Operation	59
7.2. Temperature calibration.....	62
7.3. Pumps calibration	64
7.4. Level calibration.....	66
7.5. Opening the door without supply.....	69
7.6. Failure indication (alarms)	70
8. Attachments	73
8.1. Electric scheme	73
8.2. Hydropneumatic Diagram	73
8.3. Exploded View	73
9. LogBook	73
9.1. Periodic preventive maintenance on schedule:	73
10. G.P.V. - After Sales Management	75

1. Introduction

The information contained in this manual is intended to enable users to be informed of periodic preventive maintenance in order to maintain the proper functioning of the equipment and as their execution ensures the product's limited warranty.

Also included are maintenance of safety items, a list of spare parts, calibration and maintenance routines, attachments, LogBook, among others.

Intervention by unlicensed personnel implies the termination of warranty for products still within the validity period of this warranty for the equipment or parts thereof, and THE ASSIGNMENT OF THE RESPONSIBILITY OF THE MANUFACTURER for equipment outside these periods.

Symbologies used in this manual:

Danger - Failure to follow the procedures recommended in this warning can lead to serious injury and property damage.

Attention - Failure to follow the procedures recommended in this warning may lead to equipment malfunction or property damage.

Note - The text is intended to provide important information for the correct understanding and proper functioning of the equipment.

Safety Notes:

Installation must be performed by qualified and qualified personnel.

For any reason, any safety components of the equipment must be modified or changed.

Installation must be performed by qualified and qualified personnel. Improper use of this function may impair the washing and disinfection process or cause injury to the staff.

When the equipment is hot, care must be taken when handling all unprotected internal parts of the unit, as these parts can cause burns during maintenance.

According to ANVISA regulations (see regulations specific to the country of destination) the test must be performed at the end of all preventive and corrective maintenance.

2. Baumer Limited Warranty Term

I. Introduction

This term sets forth the limited warranty conditions of the BAUMER Product to the original Purchaser, in which Baumer S.A. warrants that its products are delivered free of manufacturing defects.

II. Warranty Period

The warranty starts from the issuance of the Sales Invoice and has the following deadlines:

- 5 (five) years for pressure vessels (steam sterilizers and steam generation boilers);
- 13 (thirteen) months against manufacturing defect for general parts, not indicated in specific warranty;
- 6 (six) months for electrical / electronic material, hydraulic seals and other items listed in the specific section of the instruction manual or equipment manual;
- 90 (ninety) days for components susceptible to normal wear and tear such as lamps, plugs, fittings, connectors, instrumentation cables, removable accessories, reverse osmosis membranes and other items listed in the specific section of the instruction manual or use of the equipment;

III. Condition for granting the Limited Warranty

It is a condition for granting the warranty of the product the formal proof of compliance with the following criteria:

1. The product has been received by the Buyer, attesting to your conference with the one informed on the invoice at the time of delivery, free of apparent damages and with the packaging intact.
2. Installation has occurred as determined in the Product Installation Guide or the Installation Manual, with the correct fulfillment and approval of the particular Technical Delivery Report (RET).
3. The installation has been done by a professional accredited or formally authorized by Baumer.
4. The Installation Qualification (QI) and Operational Qualification (QO) protocols have been implemented and approved as required by law and / or mandated in the Installation Guide or Manuals / Instruction for Use.
5. Check the preventive maintenance, indicated in the Maintenance Manual, through registers;
6. Use original parts and / or components for the Equipment, that is, only those supplied by Baumer S.A.;
7. Do not allow interventions by unauthorized technical agents for repairs, applications and installations of additional components;
8. To return formally to Baumer S.A. (Mogi Mirim factory), through its authorized agent, within a maximum of 10 (ten) working days, the components and / or parts replaced as warranty.

Note: The guarantee may be suspended at any time in case of pending clauses III item 8. If there are open return parts the corrective maintenance service will be done under collection, with presentation of linked budget.

IV - Effectiveness of the Warranty

The granting of the Warranty will only be formally given to the Buyer with proof of compliance with the items in section III of this document and must be sent to Baumer S.A. at the physical or electronic address included in the Installation Guide and / or Instructions for Use or Product Manual the following documents:

1. Condition for granting the Limited Warranty

The. Receipt of the product by the shipping company with the acceptance and checked by the recipient;

B. Technical delivery report completed by the Baumer S.A. representative or by a person duly authorized for its execution;

c. Completed and approved Installation Qualification and Operation Qualification Protocols (where applicable).

2. Condition for granting the Limited Warranty

a. The proof of the performance of preventive maintenance in accordance with the Maintenance Manual.

Note: the warranty of the exchanged components is conditioned to the term of Section III.8. Upon receipt of the replaced component under warranty and, after technical analysis, if undue replacement is proven (for misuse, non-original component or situation that does not characterize product failure or warranty component), Baumer S.A. reserves the right to charge the replaced component to the Buyer or person responsible for the exchange.

V. Exclusions

The following items, among others compatible with the above, are not covered by the warranty:

1. Spot installation services;
2. Consulting services;
3. Components external to the Product;
4. Materials for cleaning, conservation and wear and tear from normal use;
5. Paper and glass;
6. Preventive maintenance workforce;
7. Qualification and process validation tests;
8. Periodic gauging and calibration of measuring and control instruments;
9. Software update of the controller (when applicable), except in cases where the program's proven failures adversely affect the operating and safety conditions;
10. Travel expenses and technician stays, freight, packaging and insurance;
11. Costs with outsourcing of processes due to corrective and preventive maintenance;
12. Damage caused by faults in water supplies, electricity (interruption, sub or over voltage, transients) or other changes subsequent to the installation qualification;
13. Damage caused by misuse, abuse, fall, neglect, recklessness or malpractice;
14. Damage caused by storage or use under conditions other than the specifications contained in the Instructions for Use or Product Manual;
15. Damage caused by equipment which produces or induces electromagnetic interference outside the regulatory limits;
16. Damage caused by accessories and Third Party Products added to a Product marketed by Baumer S.A.;
17. Damage caused by violation of the Product, attempted repair or adjustment by third parties not authorized by Baumer S.A.;
18. Damage caused by agents of nature, such as electric discharges (lightning), floods, fires, landslides, earthquakes, etc;
19. Losses and damages caused by the Product or product performance, including but not limited to loss of profits, financial losses and productivity limitations, resulting from acts related to the hypothesis of not covering this warranty;
20. Damage caused to the installed product after the expiration of the warranty periods described above;
21. Damage caused to the Product due to transportation.

VI. Manufacturer's Limitations of Liability

1. The obligations assumed by Baumer S.A. as a result of this Warranty Term are limited to those expressly included herein;
2. The solutions provided in this Warranty Term are the only and exclusive remedies offered to the customer;
3. Under no circumstances shall Baumer S.A. be liable for any direct, indirect, including loss of profit, special, incidental or consequential damages, whether based on contract, legal action, damages or other legal theory;
4. In no event shall Baumer S.A.'s liability for material damages exceed the maximum price limit of the Product that caused such damage.

VII. Additional Warranties

The extended and / or special warranties will be negotiated separately between Baumer S.A. and the Buyer. After hiring, they will be registered in a specific supply contract for this purpose.

VIII. Others

In accordance with the policy of continuous development, Baumer S.A. reserves the right to make, without prior notification, changes to the product mentioned in this document.

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without the prior permission of Baumer S.A..

3. Replacement of Safety Items

This section describes the procedures for replacing the equipment safety items.

3.1. Replacing the HMI battery

The battery exchange is informed by the control HMI screen when it is low and must be replaced. The battery is on the back of the HMI.

To replace the battery proceed as follows:

- With the equipment on (not running cycles) open the front trim using the key supplied with the equipment;
- Open the battery compartment on the back of the HMI in the upper left and remove the discharged battery;
- Insert the new model CR-2032 battery and close the compartment;
- Close the front trim.

Figure 26 - Back view of the HMI

Figure 27 - Battery compartment open

Replacing the battery with the equipment switched off will result in the loss of information recorded on the HMI.

Battery replacement should be performed by trained personnel, because the equipment is energized.

3.2. Replacing the fuses

The fuses are of glass with current of 0.5 A, they have the function of protect the controller.

To replace the fuses proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Pull the electric slide control box out until the end of the stroke, unlock and open the door;
- Locate the faulty fuse and remove it;
- Insert a new fuse and return to the above steps to finish.


Figure 28 - Replacing fuses

To identify which fuse to replace, check the wiring diagram.

3.3. Replacing the emergency button

The emergency button is located on the front of the loading side and on the two front loading side and the unloading side (when two-door equipment). Once triggered, the emergency button will stop powering the equipment.

To replace the emergency button proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Access the back of the equipment panel and remove the connection from the emergency button wires;
- Remove the contact block from the emergency button using a wrench, pressing on its clamp;
- Remove the locking ring positioned at the back of the emergency button by turning it counterclockwise;
- Remove the emergency button from its holder hole;
- Reposition the new emergency button on its support hole, rotate the locking ring clockwise;
- Reposition the new contact block on the button and re-connect the wiring connections;
- Turn on the equipment and test the emergency button to validate its operation.


Figure 29 - Emergency button on loading side


Figure 30 - Emergency button on unloading side

3.4. Replacing the general circuit breaker

The general circuit breaker is located on the electric control box.

Disconnect the main power before starting to change the component.

To replace the general circuit breaker proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Access to the electrical control box;
- Remove the connected cables with a philips wrench;
- Remove the general circuit breaker from the control;
- Reposition the new general circuit breaker;
- Reposition the input and output cables of the circuit breaker and redo the connections;
- Turn on the equipment and test the general circuit breaker to validate its operation.

Figure 31 - General circuit breaker

3.5. Replacing the unipolar circuit breaker

The unipolar circuit breaker is located on the electric control box.

Disconnect the main power before starting to change the component.

To replace the unipolar circuit breaker proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Pull the electric slide control box out until the end of the stroke, unlock and open the door;
- Access to the electrical control box;
- Remove the connected cables with a philips wrench;
- Using a wrench, press down the locking catch located on the bottom of the circuit breaker;
- Remove the bipolar circuit breaker from the control;
- Reposition the new circuit breaker;
- Reposition the input and output cables of the circuit breaker and redo the connections;
- Turn on the equipment and test the motor circuit breaker to validate its operation.

Figure 32 - Unipolar circuit breaker

3.6. Replacing the motor circuit breaker

The water pump motor circuit breaker is located on the electric generator control.

Disconnect the main power before starting to change the component.

To replace the motor circuit breaker proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Pull the electric slide control box out until the end of the stroke, unlock and open the door;
- Access to the electrical control box;
- Remove the connected cables with a philips wrench;
- Remove the circuit breaker by pushing it down, relieving the fastening spring;
- Reposition the new motor circuit breaker;
- Reposition the input and output cables of the circuit breaker and redo the connections;
- Turn on the equipment and test the motor circuit breaker to validate its operation.

Figure 33 - Motor circuit breaker

3.7. Replacing the air pressure switch

The air pressure switch is located at the top of the product box.

Disconnect the main power before starting to change the component.

To replace the pressure switch proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Access the air pressure switch;
- Remove the connected cables with an allen wrench;
- Remove the Air hose by pressing the quick coupler down;
- Remove the pressure switch with a wrench;
- Reposition the new pressure switch;
- Reposition the cables and air hose and redo the connections;
- Turn on the equipment and test the pressure switch to validate its operation.

Figure 34 - The air pressure switch

The replaced pressure switch must be calibrated.

3.8. Replacing the air pressure switch on the channel

The air pressure switch on the channel is located at the top of the loading and unloading doors.

Disconnect the main power before starting to change the component.

To replace the pressure switch proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Access the air pressure switch;
- Remove the connected cables with an allen wrench;
- Remove the air hose by pressing the quick coupler down;
- Remove the pressure switch with a wrench;
- Reposition the new pressure switch;
- Reposition the cables and air hose and redo the connections;
- Turn on the equipment and test the pressure switch to validate its operation.

Figure 35 - Channel air pressure switch

The replaced pressure switch must be calibrated.

3.9. Replacing the air pressure switch

The circulation pressure switch is located on the right side of the equipment.

Disconnect the main power before starting to change the component.

To replace the pressure switch proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Access the air pressure switch;
- Remove the connected cables with an allen wrench;
- Remove the pressure switch with a wrench;
- Reposition the new pressure switch;
- Reposition the cables and air hose and redo the connections;
- Turn on the equipment and test the pressure switch to validate its operation.

Figure 36 - The air pressure switch

The replaced pressure switch must be calibrated.

3.10. Replacing the anti-crushing system

The anti-crushing system is located at the top of the loading door and the option of two-door equipment also at the unload door. During the door ascending, if an obstruction occurs, the system is activated and the door returns to the open position and a visual and audible alarm is emitted at the HMI.

To replace the door micro switch proceed as follows:

- Open the equipment door by pressing the "OPEN DOOR" button on the HMI;
- Turn the equipment off and open the front trim using the key supplied with the equipment;
- Remove the two fixing screws of the door trim using an allen wrench;
- Remove the door trim;
- Disconnect the wires from the micro switch;
- Remove the fixing screws of the holder of the micro switch, with an aid of a screwdriver;
- Install the new micro switch and return the above steps to finish;
- If necessary, adjust the position of the holder and test to validate its operation.


Figure 37 - Replacing the anti-crushing system

To replace the anti-crushing system proceed as follows:

- Open the equipment door by pressing the "CLOSE DOOR" button on the HMI;
- When start the upward movement of the door, press the top of the door using the hand.

The door must stop the ascend movement and begin to descend and a visual and audible signal is emitted on the HMI.

3.11. Replacing the thermostat the water resistances

The thermostat of the water resistances is located below the loading door of the equipment.

Disconnect the main power before starting to change the component.
of the component.

To replace the thermostat proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Access the thermostat of the water resistant.
- Remove the thermostat cover with a screwdriver;
- Remove the connected cables with a screwdriver;
- Remove the fixing screws of the thermostat;
- Remove the thermostat;
- Remove the capillary of the thermostat with a screwdriver;
- Reposition the new thermostat and fasten the screws;
- Reposition the cables and the capillary of the thermostat and redo the connections;
- Close the thermostat cover;
- Turn on the equipment and test the thermostat to validate its operation.

Figure 38 - Electrical resistance thermostat

3.12. Replacing the air resistances thermostat

The thermostat of the air resistances is located below the loading door of the equipment.

Disconnect the main power before starting to change the component.
of the component.

To replace the thermostat proceed as follows:

- With the equipment off, open the front trim using the key supplied with the equipment;
- Access the air resistant thermostat.
- Remove the temperature switch cover with a screwdriver;
- Remove the connected cables with an allen wrench;

- Remove the fixing screws of the temperature switch;
- Remove the thermostat;
- Remove the capillary of the thermostat with a screwdriver;
- Reposition the new temperature switch and fasten the screws;
- Reposition the cables and the capillary of the thermostat and redo the connections;
- Close the thermostat cover;
- Turn on the equipment and test the thermostat to validate its operation.

Figure 39 - Air resistance thermostat

4. Other maintenance

4.1. Door gasket replacement and maintenance

- The life of a gasket is directly proportional to the intensity of use of the equipment (number of hours in operation, number of cycles per day).
- Never use lubricant on the gasket as it must be dry-worked by an inflatable gasket;
- To remove the gasket, never use sharp or pointed tools, such as a screwdriver, a knife, a fork, or any other. Use only the special spatula with curved tip (Baumer code 86428);
- Never remove a gasket with hot equipment, as it may rupture the splice.

4.1.1. To remove the gasket:

4.1.1.1. Disconnect the air intake from the trim located in the upper area of the chamber;

4.1.1.2. Loosen the cap nuts on the chamber insider;


Figure 40 – Gasket locks disassembly

4.1.1.3. Position the spatula between the channel and the gasket by pulling it out until it completely releases from the chamber channel and remove all the gasket;

Figure 41 – Gasket removal

4.1.2. Gasket assembly:

The gasket is provided with a membrane that serves as a backrest between the latches (Baumer code 906295) and the channel.

- 4.1.2.1. If the gasket comes with pre-defined accommodation holes for the latches in the gasket backrest membrane, skip to section **Erro! Fonte de referência não encontrada.**;
- 4.1.2.2. If the gasket does not come with the pre-defined holes, proceed as described in the section 4.1.2.3;
- 4.1.2.3. Position the gasket in its chamber channel. Proceed as described in the sections 4.1.2.8 a 4.1.2.17;
- 4.1.2.4. With a soft tip marker or pen, paint the gasket in the attachment points (chamber channel holes);
- 4.1.2.5. Remove the gasket from its channel as described in the section 0;
- 4.1.2.6. Open the gasket backrest membrane when attaching the latches, which can be done using the latches or a blade utility knife;


Figure 42 - Membrane to be opened

Figure 43 - Open Membrane

Do not puncture or damage the gasket.

- 4.1.2.7. Insert the latches on the pre-defined holes positioning or in the membrane open regions matching the chamber channel attachment positioning;


Figure 44 - Latch to be inserted


Figure 45 - Latch inserted

- 4.1.2.8. Place the air intake nozzle in the connection tube located in the upper chamber channel (position 1, as indicated in both Figure 46 and Figure 47);


Figure 46 – Air intake nozzle placement

- 4.1.2.9. Place the opposite end gasket latch opposite in the lower chamber channel central position (position 2, as indicated in the Figure 47);
- 4.1.2.10. Stretch the gasket from the upper and lower attachment distributing the gasket uniformly in a way that it stay placed in the chamber channel bottom surface and locked with its latches (position 3, as indicated in the Figure 47);


Figure 47 – Gasket placement sequence

- 4.1.2.11. Screw the gasket guide (Baumer code 906551) in the gasket latch;
- 4.1.2.12. Insert the guide into the hole and house the gasket on the channel;


Figure 48 - Gasket lock positioning on the channel hole

- 4.1.2.13. Release the guide;
- 4.1.2.14. Fix the latches through the cap nuts;
- 4.1.2.15. Repeat the steps as described in the sections from 4.1.2.11 until 4.1.2.14 for each latch;
- 4.1.2.16. After the gasket assembly, check if there is not any remnant parts outside the chamber channel, especially in the upper and lower part of the channel;

Figura 49 - Verify any remnant parts of the gasket in the chamber channel

- 4.1.2.17. Connect the compressed air intake to the gasket;

The air pressure applied to the gaskets shall be between 0,8 bar and 1,0 bar;

- 4.1.2.18. Switch on the equipment and test the gasket operation (check for signs of air leakage in the gasket or water leakage in the door sealing system).

- 4.1.2.19.

The gap between the door and the equipment chamber shall be between 3 and 5 at any region.
For further details regarding the door adjustment contact Baumer After Sales Management.

5. Preventive Maintenance Schedule

5.1. Periodic preventive maintenance on schedule:

Failure to do so will result in loss of warranty.

5.1.1. Quarterly:

- Examine the filter of the water inlet manifold, if there is a change in the flow rate or the initial color of the filter, proceed with the exchange;
 - Examine all piping for leaks. If leaks are found, perform the re-tighten. If leakage continues, then disassemble, apply tape seal PTFE thread and re-tighten. If the leak persists, then perform the item exchange;
 - Retight the electrical power contacts;
 - Check the air supply pressure of the gaskets, which can vary from 0.8 to 1.0 bar;
 - Check the adjustment between the door (s) and the channel (s). If it is outside the range from 3 to 5mm at any region make the adjustment;
 - Inspect the doors gaskets, if any faults (cuts, deformations or other) are found to proceed with the exchange;
 - Carry out the cleaning of the resistances of the chamber, boiler and fan, using a brush of stainless steel;
 - Clean Y filters.
 - Clean the showers by disassembling them, if there is encrusted material, spray on the internal and external surface the FACTS® descaling (Baumer HER-501 Code).
 - Check all hydraulic hoses for leaks. If leaks are found, perform the hose exchange. If leakage is at the hose connection, re-tighten the connection and if necessary replace them;
-

When performing a repair on hoses or pipes, be aware of the presence of liquids inside the hoses or pipes.

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5.1.2. Semiannually:

- Replace the aeration filter (or every 1,500 cycles, whichever comes first);
- Replace door gasket;
- Clean the seat of the check valves;
- Replace the water inlet manifold filter.

5.1.3. Annually:

- Calibrate all temperature sensors, pressure switches and pressure sensors;
- Perform the Re-qualification of Installation, Operation and Performance.

6. Spare parts

DESCRIPTION	BAUMER CODE
Silicone ring 2"	55301
42.5 x 2.6 O'ring	891228
Side injector nozzle E-2002	890471
Auxiliary circuit-breaker motor	891848
Vertical level	892009
Pump 4.0cv in brass 50Hz - 1b	901277
Pump 4.0cv in brass 60Hz - 1b	890849
Single-stage centrifugal pump 3 hp 50Hz - 2b	901275
Single-stage centrifugal pump 3 hp 60Hz - 2b	904054
Chemical nickel pump 0.5hp 50 Hz - 1b	901278
Chemical nickel pump 0.5hp 60 Hz - 1b	893902
Thermo Supply Water Pump 50 Hz - 2b	901276
Thermo Supply Water Pump 60Hz - 2b	904055
Peristaltic pump TW	903540
Emergency button push type	98278
White light button	99013
Bell	905843
Stainless steel level switch	903641
On / Off switch	88256
Pneumatic cylinder 32 x 720 - Term 287	901494
Pneumatic cylinder 32 x 780 - Term 400	895636
Set Check valve 1/2"	86154
25A 24VCC Contactor	55047
8DI + 6DO PLC relay controller	901554
Unipolar circuit breaker 10A	88459
6 to 10A Motor circuit breaker	891466
2.5 to 14A Motor circuit breaker	891196
Tripolar circuit breaker 40A	88458
Tripolar circuit breaker 63A	88456
4DI PT-100 Temperature Measurement Module	901559
8DI + 6DO Digital Expansion Module	901555
Filter "Y" 1/2 "	45062
HEPA Absolute filter 305x305x78	891053
Switching Power Supply 24V, 6.25A	902220
Custom Printer Power Supply	89548
Boiler gasket TW-400	906256
Boiler gasket TW-287	906255
Boiler bottom gasket	906254
Inflatable gasket TW-287	906266
Inflatable gasket TW-400	906267


DESCRIPTION	BAUMER CODE
TW lamp gasket	903601
Boiler upper gasket	906253
Touch screen HMI 5.7"	901556
Custom serial printer	892739
Plus-s4b Thermal Printer	896504
Led dichroic MR16 lamp	903115
Hydraulic Water Hose ¾ "300 psi	101097
Hydraulic connection hose with pump 2" 300 psi	892554
Micro Fan 220VAC	87369
Miniature Switch (Leaf spring lever with roller)	900361
Mini Pressure Regulator	89974
Spring inner diameter 31 x 43 length	890547
Pre filter EU4 305x305x15	891054
Pressure switch 0.1 to 1.0 bar NO	88442
Pressure switch 3 levels E-2002	891100
Pressure switch NO	88441
Pressure switch NO	88442
Pressure regulator 0 to 3.5 BAR	898728
16A 24VCC relay - 1 reversible contact	905010
Resistance 220V 5000W	906265
Boiler Resistance 4000W 220V	905515
Bottom resistance 900W 220V	892025
Bigger resistance 1200W 220V	891160
Medium Resistance 1200W 220V	891159
Lower resistance 1200W 220V	891158
Inductive Sensor M12 PNP NA	899161
Temperature sensor type PT-100	895714
Temperature sensor type PT-100 (double)	905695
MR16 Led Socket bulb	903116
Safety thermostat 170°C	900984
Butterfly valve 2 clamp	890948
Solenoid valve ½" water - 2/2 NC brass	897412
Solenoid valve 3/8" water - 2/2 NC stainless steel	894416
Retention valve - ½" NPT brass	55375
Non-return valve – 6mm O/D Tube	892401
Fan Siroco thermo 50Hz	900255
Fan Siroco thermo 60Hz	891056
Round glass lamp TW	903117

7. Calibration and Maintenance Routines


The washer-disinfector of TW-E2000 line has functions for maintenance procedures by the maintenance team. The main functions are described below:

7.1. Manual Operation

The product is designed to operate fully automatically. To assist maintenance and check the correct operation of the drives, the product can be operated manually, through the access password specified by the supervisor. (Note: User registration information and access password are contained in the Instructions for Use).


In the start menu press the "CONFIG" key.

The screen with "SUPERVISOR", "PASSWORD", "FACTORY", "LOGOFF", "RETURN", will open.

Pressing the field below the "FACTORY" text box will display an alphanumeric keypad, enter the default access password and confirm.


Click on the option below the password edit field to access the "FACTORY" screen.

The screen with the options "OPTIONAL", "ACCESS RESTRICTED", "LOAD DEFAULT", "OPERATION MANUAL", "GRAPHICS", and "RETURN" will open.


Pressing the field below the "OPERATION MANUAL" text box will display an alphanumeric keypad, enter the default access password and confirm.

Click on the option below the password edit field to access the "OPERATION MANUAL" screen.

The screen with the "SYSTEM WATER", "SYSTEM PRODUCT", "VARIOUS", and "SYSTEM DOORS" options will open.


Press the "WATER SYSTEM" text box.


The "WATER SYSTEM" screen allows access to the control buttons, "COLD WATER VALVE", "HOT WATER VALVE", "TREATED WATER VALVE", "WATER PUMP" and "DRAINAGE VALVE", "CHAMBER WATER VALVE", and / or "BOILER WATER VALVE", to make the accesses as needed.

Example: To check the water levels, press "DRAINAGE VALVE", "COLD WATER VALVE" and wait until the levels are reached.

To return to the "MANUAL OPERATION" screen, press the icon that symbolizes an "OPEN DOOR".

In the "MANUAL OPERATION" screen, press the "SYSTEM PRODUCT" text box to access the product entry through the peristaltic pumps.


The "PRODUCT SYSTEM" screen allows access to the control buttons, "PRODUCT 1" and "PRODUCT 2", "PRODUCT 3", and "PRODUCT 4".


Note: The number of products enabled depends on the number of pumps purchased at the time of purchase.

Example: If you have purchased 2 peristaltic pumps, "PRODUCT 1" and "PRODUCT 2" will be enabled and the remaining ones will be visible, but disabled.

To return to the "MANUAL OPERATION" screen, press the icon that symbolizes an "OPEN DOOR".


In the "MANUAL OPERATION" screen, press the "VARIOUS" text box to access the Exhaust Fan, Water Circulation Pump, Audible Alarm, Camera Illumination and End of Cycle Bright Button.


The "VARIOUS" screen allows on / off access to the "AIR CIRCULATION" and "WATER CIRCULATION", "CIRCULATION VALVE", "INTERNAL LIGHT", "UNLOAD DOOR LIGHT" and "BUZZER" access as needed.

Attention: To preserve the mechanical seal life of the "WATER CIRCULATION" Pump, even if the Pump is activated, it will only be switched on if one of the levels is reached.

To return to the "MANUAL OPERATION" screen, press the icon that symbolizes an "OPEN DOOR".


In the "MANUAL OPERATION" screen, press the "SYSTEM DOORS" text box to access the doors controls.


The "SYSTEM OF DOORS" screen allows access for opening / closing the doors "CLOSE LOADING DOOR" and "CLOSE UNLOADING DOOR".

Attention: To preserve the interlock system between the ports, the software does not allow the openings of both ports simultaneously.

To return to the "MANUAL OPERATION" screen, press the icon that symbolizes an "OPEN DOOR".


In the "MANUAL OPERATION" screen, press the "RETURN" text box to access the previous menu.

The screen with the options "OPTIONAL", "ACCESS RESTRICTED", "LOAD DEFAULT", "OPERATION MANUAL", "GRAPHICS", and "RETURN" will open.


Press the "RETURN" text box to access the previous menu.

The screen with "SUPERVISOR", "PASSWORD", "FACTORY", "LOGOFF", "RETURN", will open.

Press the "RETURN" text box to access the "START" menu.

7.2. Temperature calibration


In the start menu press the "CONFIG" key.


The screen with "SUPERVISOR", "PASSWORD", "FACTORY", "LOGOFF", "RETURN", will open.

Press the blank field below the "SUPERVISOR" text box.

The screen with "SUPERVISOR", "PASSWORD", "FACTORY", "LOGOFF", "RETURN", will open.

Pressing the field below the "CALIBRATION" text box will bring up the access control screen.


Pressing the blank field in front of the "User Name" text box will pop up an alphanumeric keypad,

Enter the user name using the alphanumeric keypad and confirm by pressing the "ENT" key;

Pressing the blank field in front of the "Password" text box will bring up an alphanumeric keypad;


Enter the Administrator password using the alphanumeric keypad and confirm by pressing the "ENT" key.
The "SENSOR CALIBRATION" screen will appear.

The screen with "SUPERVISOR", "PASSWORD", "FACTORY", "LOGOFF", "RETURN", will open.

Pressing the field below the "CALIBRATION" text box will bring up the access control screen.

The "CALIBRATION" option allows the calibration of the TE-6A, TE-6B and TE-9 sensors.


To calibrate the temperature sensors, it will be necessary to use a calibration oven set at a reference temperature for comparison. Proceed as below to calibrate the temperature sensors:


- Turn on the Calibration Oven and insert the Sensor to be calibrated.
- Set the Calibration Oven temperature to 80°C ± 0.5°C.
- Check the reading of the sensor to be calibrated with the calibration oven indication.
- If there is a difference in temperature, adjust as below.
- Calculate the difference between the reading of the Sensor to be calibrated and the Calibration oven.
- The value found must be increased if the result is positive or decremented, if negative, of the standard factor of adjustment. The default value of "AJ" is 15.
- To adjust, press on the parameter to be modified, enter the desired value and confirm with the "ENTER" key.
- Check again the reading of the control with the indication of the oven and if necessary adjust again.

To exit press the "RETURN" buttons.

7.3. Pumps calibration


In the start menu press the "CONFIG" key.

The screen with "SUPERVISOR", "PASSWORD", "FACTORY", "LOGOFF", "RETURN", will open.

Press the blank field below the "SUPERVISOR" text box.

The screen with "SUPERVISOR", "PASSWORD", "FACTORY", "LOGOFF", "RETURN", will open.

Pressing the field below the "CALIBRATION" text box will bring up the access control screen.

Pressing the blank field in front of the "User Name" text box will pop up an alphanumeric keypad,


Enter the user name using the alphanumeric keypad and confirm by pressing the "ENT" key;

Pressing the blank field in front of the "Password" text box will bring up an alphanumeric keypad;

Enter the Administrator password using the alphanumeric keypad and confirm by pressing the "ENT" key.


The "SENSOR CALIBRATION" screen will appear.

The screen with "CALIBRATION", "PUMPS", and "LEVEL" will open.

Pressing the field below the "PUMPS" text box will bring up the access pumps calibration.

The "PUMPS" option allows calibration of the peristaltic pumps that control the intake of product.


Calibration should be performed on all peristaltic pumps added to the equipment, and may be from one (01) to four (04) pumps. To perform the calibration of the peristaltic pumps will require the use of a Beaker of 1000ml complete with water. Proceed as below to calibrate the peristaltic pumps:


- Insert the pump suction device into the Becker containing 1000ml.
- Press the "NEW" key.
- Select the pump to be calibrated from 01 to 04.
- In the "DEFINE VOLUME USED" text box, enter 100ml.
- Press the "START" key, and check the consumption of 100ml in the Becker and press the "STOP" key.
- Press the "SAVE" key.
- Repeat the procedure until all available pumps are calibrated.

To exit press the "RETURN" buttons.

7.4. Level calibration


In the start menu press the "CONFIG" key.


The screen with "SUPERVISOR", "PASSWORD", "FACTORY", "LOGOFF", "RETURN", will open.

Press the blank field below the "SUPERVISOR" text box.

The screen with "CALIBRATION", "DATE AND TIME", "INSTITUTE NAME", and "RETURN" will open.

Pressing the field below the "CALIBRATION" text box will bring up the access control screen.

Pressing the blank field in front of the "User Name" text box will pop up an alphanumeric keypad,


Enter the user name using the alphanumeric keypad and confirm by pressing the "ENT" key;

Pressing the blank field in front of the "Password" text box will bring up an alphanumeric keypad;


Enter the Administrator password using the alphanumeric keypad and confirm by pressing the "ENT" key.

The "SENSOR CALIBRATION" screen will appear.

The screen with "CALIBRATION", "PUMPS", and "LEVEL" will open.

Pressing the field below the "CALIBRATION" text box will bring up the access control screen.


The "LEVEL PRESSURE SWITCH" option allows the calibration of the "LOW LEVEL", "MEDIUM LEVEL" and "HIGH LEVEL".

Calibration shall be performed for all three water levels using a calibrated flow meter, so that the low, medium and high level shall meet the requirements of the table below:

Level	TW-E2000-287(P)	TW-E2000-400(P)
Low (l)	25	28
Medium (l)	30	35
High (l)	45	48

To carry out the calibration of the pressure switch levels, proceed as follows:

Low Level Calibration:

- Select the "LOW LEVEL" text box.
- Loosen the three pressure level adjustment screws.
- Press the button "COLD WATER VALVE".
- The water inlet is released until it reaches the predetermined level.
- If the specified level is not as specified, press the "OPEN DRAIN VALVE" button, so that the chamber is completely emptied.
- Adjust the adjusting screw in order to get more or less water, depending on the previous result.
- Repeat the previous steps until the volume of the specified level conforms to the table.
- Perform the procedure 5 times to ensure repeatability of the process.

Medium Level Calibration:

- Select the "MEDIUM LEVEL" text box.
- Press the button "COLD WATER VALVE".
- The water inlet is released until it reaches the predetermined level.
- If the specified level is not as specified, press the "OPEN DRAIN VALVE" button, so that the chamber is completely emptied.
- Adjust the adjusting screw in order to get more or less water, depending on the previous result.
- Repeat the previous steps until the volume of the specified level conforms to the table.
- Perform the procedure 5 times to ensure repeatability of the process.

High Level Calibration:

- Select the "HIGH LEVEL" text box.
- Press the button "COLD WATER VALVE".
- The water inlet is released until it reaches the predetermined level.
- If the specified level is not as specified, press the "OPEN DRAIN VALVE" button, so that the chamber is completely emptied.
- Adjust the adjusting screw in order to get more or less water, depending on the previous result.
- Repeat the previous steps until the volume of the specified level conforms to the table.
- Perform the procedure 5 times to ensure repeatability of the process.


The screens are merely illustrative with variations due to the options configurations requested at the time of purchase.

7.5. Opening the door without supply


This operation should only be carried out by trained technicians.

In the need to open the doors in the absence of supplies, compressed air and electricity, follow the instructions below:

Open load door:

- Open the frontal load trim and access the valve block;

- On the SYV-5 valve using a screwdriver, position the upper manual override button by pressing it and turning it to a vertical position;

- The channel will be depressurized.

- On the SYV-3 valve and using a screwdriver, position the top manual override button by pressing it and turning it to a vertical position;

- The door cylinders will be depressurized and the door will be open.

To return for initial position:

- Reposition the SYV-5 valve manually in the horizontal position;
Reposition the SYV-3 valve manually in the horizontal position.

To open the unload door:

- Follow the procedures above changing the valve SYV-3 to SYV-4.

If the drive buttons are not positioned in the original factory position, the machine will not operate.

7.6. Failure indication (alarms)

TW-E2000 washer-disinfector line have an extremely functional control logic so that they do not allow actions or parameters outside the required ranges or conditions, which makes this product extremely reliable during the washing disinfection cycle.

The washer-disinfector alarms may occur under different circumstances or process steps.

The alarms are different from the other parameters because they are flashing intermittently on the screen.


Figure 50 - Failure indication

See below the main alarms and their causes and solutions:

Problem	Cause	Corrective measure
High water level failure.	• Level pressure switch.	• Evaluate level pressure switch. • Adjust the average pressure switch level.
Supply Failure.	• Cold water inlet. • Hot water inlet. • Purified water inlet. • Compressed air inlet. • Steam inlet (when steam heating).	• Check if there is cold water and if it is connected to the point of the equipment. • Check if there is hot water and if it is connected to the point of the equipment. • Check for purified water and if it is connected to the point of the equipment. • Check if there is compressed air and if it is connected to the point of the equipment. • Check if there is steam and if it is connected to the point of the equipment.
Pump circulation failure.	• Circulation pump	• Check motor circuit breaker. • Check that the pump is mechanically locked. • Check that the pump motor is not locked.
Temperature sensor failure.	• Thermo resistance (PT100).	• Check or replace the temperature sensor.
Cold water faulty.	• Cold water inlet.	• Check if there is cold water and if it is connected to the point of the equipment.
Hot water faulty.	• Hot water inlet.	• Check if there is hot water and if it is connected to the point of the equipment.
Purified water faulty.	• Purified water inlet.	• Check for purified water and if it is connected to the point of the equipment.
Lack of Compressed air.	• Compressed air inlet.	• Check if there is compressed air and if it is connected to the point of the equipment.
Lack of steam.	• Steam inlet.	• Check if there is steam and if it is connected to the point of the equipment.
Close the door and then start the cycle.	• Loading door	• Close the load door and then start the cycle. • Check inductive sensor that indicates closed door. • Check that there is an obstacle to closing the door.

Problem	Cause	Corrective measure
Close unload door.	• Unload door	• Close the unload door to release the load door. • Check inductive sensor that indicates closed door. • Check that there is an obstacle to closing the door.
Product 01 has reached the level of safeness. Refilling.	• Container with Product 01. • Float level.	• Check product container 01 and refill if necessary. • Check that the float magnetic sensor is working.
Product 02 has reached the level of safeness. Refilling.	• Container with Product 02. • Float level.	• Check product container 02 and refill if necessary. • Check that the float magnetic sensor is working.
Product 03 has reached the level of safeness. Refilling.	• Container with Product 03. • Float level.	• Check product container 03 and refill if necessary. • Check that the float magnetic sensor is working.
Product 04 has reached the level of safeness. Refilling.	• Container with Product 04. • Float level.	• Check product container 04 and refill if necessary. • Check that the float magnetic sensor is working.
Failure of the heating air system	• Safety thermostat of heating air system. • Drive system (controller output, contactor, relay and fuse).	• Check that the thermostat is not open. • Check for any drive component failure.
Failure of the heating water system	• Safety thermostat of heating water system. • Drive system (controller output, contactor, relay and fuse).	• Check that the thermostat is not open. • Check for any drive component failure.
Failure of compressed air in the loading door channel.	• Pressure drop in the loading door sealing system.	• Check that the inflatable gasket is not punctured. • Check that the packing is connected to the pressure point, • Check that the supply pressure of the gasket is between 0.8 and 1.0 bar.
Failure of compressed air in the unloading door channel.	• Pressure drop in the unloading door sealing system.	• Check that the inflatable gasket is not punctured. • Check that the packing is connected to the pressure point, • Check that the supply pressure of the gasket is between 0.8 and 1.0 bar.
Equipment operating without Rack (only when equipment operating with two pumps).	• Cycle started without rack (central showers without operation).	• Abort the cycle and load the rack and then restart the cycle.

8. Attachments

8.1. Electric scheme

Attached at the end of this manual, code:
For Washer-disinfector without boiler or with two pumps: E-TW-1.
For Washer-disinfector with boiler: E-TW-2.

8.2. Hydropneumatic Diagram

Attached at the end of this manual, code:
For Washer-disinfector without boiler or with two pumps: H-TW-01-p.
For Washer-disinfector with boiler: H-TW-02-p.

8.3. Exploded View

Attached at the end of this manual, code: V-TW-01-p.

9. LogBook

The execution of all periodic preventive maintenance presented above is the essential condition for maintaining the limited warranty according **Baumer Limited Warranty Term**.

Below are the fields for recording all periodic preventive maintenance performed on the equipment.

 Require that all periodic preventive maintenance be recorded in their respective fields.

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Exemplo:

Y3 - SY3 (183) → N° do cabo
 → TAG do componente
 → Entrada ou saída

A	B	C	D	E	F	G	H	I	J
Thermal Washer Disinfector									
Digital inputs					Saídas digitais				
CPU: DVP14SS211R X.0 - SM-3A (100) - Load door sensor; X.1 - SM-4A (101) - Unload door sensor; X.2 - SM-3B e SM-3C (102) - Load door safety device; X.3 - SM-4B e SM-4C (103) - Unload door safety device; X.4 - PS-5 (104) - Doors seal pressure switch; X.5 - Reserved; X.6 - S3 (106) - Unload door button; X.7 - Reserved.					CPU: DVP14SS211R Y.0 - SYV-5A (50) - Door seal pressure release valve; Y.1 - SYV-5B (51) - Door seal air valve; Y.2 - SYV-3A (52) - Opened load door valve; Y.3 - SYV-3B (53) - Closed load door valve; Y.4 - SYV-4A (54) - Opened unload door valve; Y.5 - SYV-4B (55) - Closed unload door valve.				
Expansão 1: DVP16SP11R X2.0 - PS-6AL (107) - Chamber low level; X2.1 - PS-6AM (108) - Chamber middle level; X2.2 - PS-6AH (109) - Chamber high level; X2.3 - LSRH-8A (110) - Low level product 1; X2.4 - LSRH-8B (111) - Low level product 2; X2.5 - LSRH-8C (112) - Low level product 3; X2.6 - LSRH-8D (113) - Low level product 4.					Expansão 1: DVP16SP11R Y2.0 - SY-6A (56) - Drain chamber valve; Y2.1 - SYV-2C (57) - Cold water inlet valve; Y2.2 - SYV-2B (58) - Hot water inlet valve; Y2.3 - SYV-2A (59) - Treated water inlet valve; Y2.4 - P-8A (60) - Product 1 dosing pump; Y2.5 - P-8B (61) - Product 2 dosing pump; Y2.6 - P-8C (62) - Product 3 dosing pump; Y2.7 - P-8D (63) - Product 4 dosing pump.				
Expansão 2: DVP16SP11R X3.0 - Q3 (114) - Auxiliary contact circulation pump circuit breaker; X3.1 - PS-1 (115) - Air pressure switch; X3.2 - PS-6C (116) - Low pressure indication pressure switch; X3.3 - PS-6B (117) - High pressure indication pressure switch; X3.4 - Q4 (118) - Auxiliary contact circulation and supply pump circuit breaker.					Expansão 2: DVP16SP11R Y3.0 - K1 (84) - Water heating contactor; Y3.1 - K2 (85) - Air heating contactor; Y3.2 - K3 (66) - Circulation pump contactor; Y3.3 - RL1 (67) - Water feed pump relay; Y3.4 - RL2 (68) - Drying relay; Y3.5 - RL3 (69) - Lamp in chamber relay; Y3.6 - H2 (70) - Unload door lamp; Y3.7 - H3 (71) - Buzzer.				
Expansão 3: DVP16SP11R X4.0 - SF-8 (122) - Product flow sensor.					Expansão 3: DVP16SP11R Y4.0 - SY-6C (72) - Circulation water valve (optional); Y4.1 - SYV-9B (73) - Water inlet valve to chamber (optional).				
Entradas analógicas									
Expansão 1: DVP04PT CH1 - TE-6A - Control temperature sensor.									
Expansão 2: DVP04PT CH1 - TE-6B - Verification temperature sensor.									

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


Exemplo:


Y3 - SY3 (183) → N° do cabo
 ↳ TAG do componente
 ↳ Entrada ou saída

A	B	C	D	E	F	G	H	I	J
<u>Thermal Washer Disinfector</u>									
Digital inputs					Saídas digitais				
CPU: DVP14SS211R X.0 - SM-3A (100) - Load door sensor; X.1 - SM-4A (101) - Unload door sensor; X.2 - SM-3B e SM-3C (102) - Load door safety device; X.3 - SM-4B e SM-4C (103) - Unload door safety device; X.4 - PS-5 (104) - Doors seal pressure switch; X.5 - Reserved; X.6 - S3 (106) - Unload door button; X.7 - Reserved.					CPU: DVP14SS211R Y.0 - SYV-5A (50) - Door seal pressure release valve; Y.1 - SYV-5B (51) - Door seal air valve; Y.2 - SYV-3A (52) - Opened load door valve; Y.3 - SYV-3B (53) - Closed load door valve; Y.4 - SYV-4A (54) - Opened unload door valve; Y.5 - SYV-4B (55) - Closed unload door valve.				
Expansão 1: DVP16SP11R X2.0 - PS-6AL (107) - Chamber low level; X2.1 - PS-6AM (108) - Chamber middle level; X2.2 - PS-6AH (109) - Chamber high level; X2.3 - LSRH-8A (110) - Low level product 1; X2.4 - LSRH-8B (111) - Low level product 2; X2.5 - LSRH-8C (112) - Low level product 3; X2.6 - LSRH-8D (113) - Low level product 4.					Expansão 1: DVP16SP11R Y2.0 - SY-6A (56) - Drain chamber valve; Y2.1 - SYV-2C (57) - Cold water inlet valve; Y2.2 - SYV-2B (58) - Hot water inlet valve; Y2.3 - SYV-2A (59) - Treated water inlet valve; Y2.4 - P-8A (60) - Product 1 dosing pump; Y2.5 - P-8B (61) - Product 2 dosing pump; Y2.6 - P-8C (62) - Product 3 dosing pump; Y2.7 - P-8D (63) - Product 4 dosing pump.				
Expansão 2: DVP16SP11R X3.0 - Q3 (114) - Auxiliary contact circulation pump circuit breaker; X3.1 - PS-1 (115) - Air pressure switch; X3.2 - PS-6C (116) - Low pressure indication pressure switch; X3.3 - PS-6B (117) - High pressure indication pressure switch; X3.4 - Reserved; X3.5 - Reserved; X3.6 - LSRS-9 (120) - High level boiler; X3.7 - LSRH-9 (121) - Low level boiler.					Expansão 2: DVP16SP11R Y3.0 - K1 (84) - Water heating contactor; Y3.1 - K2 (85) - Air heating contactor; Y3.2 - K3 (66) - Circulation pump contactor; Y3.3 - RL1 (67) - Water feed pump relay; Y3.4 - RL2 (68) - Drying relay; Y3.5 - RL3 (69) - Lamp in chamber relay; Y3.6 - H2 (70) - Unload door lamp; Y3.7 - H3 (71) - Buzzer.				
Expansão 3: DVP16SP11R X4.0 - SF-8 (122) - Product flow sensor.					Expansão 3: DVP16SP11R Y4.0 - Reserved; Y4.1 - SYV-9B (73) - Water inlet valve to chamber; Y4.2 - K5 (74) - Boiler contactor; Y4.1 - SYV-9A (75) - Water inlet valve to boiler; Y4.1 - SY-9 (76) - Boiler discharge boiler.				
Entradas analógicas									
Expansão 1: DVP04PT CH1 - TE-6A - Control temperature sensor; CH2 - TE-9 - Boiler temperature sensor.									
Expansão 2: DVP04PT CH1 - TE-6B - Verification temperature sensor.									
					To help understand: Example: <pre> graph LR Tag[Y3 - SY3 (183)] --> Wire[Wire number] Tag --> Component[Component TAG] Tag --> IO[Input or output] </pre>				

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